

# MOHS Kryostat

- Each cut about 8  $\mu\text{m}$
- Temperature about -23 degrees
- Distance between each cut 150  $\mu\text{m}$
- Fixation for about 10 sec in Delaunay  
(250 ml Acetone + 250 ml Alkohol 100% + 150  $\mu\text{l}$  Trichloracetic)
- Drying
- Hämalauun 1 min
- Washing in about 5 portions of water
- Eosin 40 sec
- Washing in water
- Washing in about 4 portions of alcohol 100%
- Washing in Xylol
- Stock up with a glas and Eukitt